COMPUTED TOMOGRAPHY STUDY: SIMULTANEOUS FAT AND MUSCLE EFFECT

COMPUTED TOMOGRAPHY (CT) BASED EVIDENCE OF SIMULTANEOUS CHANGES IN HUMAN ADIPOSE AND MUSCLE TISSUES FOLLOWING A HIGH INTENSITY FOCUSED ELELCTRO-MAGNETIC FIELD (HIFEM®) APPLICATION: A NEW METHOD FOR NON-INVASIVE BODY SCULPTING.


David E. Kent M.D.¹, Carolyn I. Jacob M.D.²

1. Dermatologic Surgery Specialists, Macon GA, USA; 2. Chicago Cosmetic Surgery and Dermatology, Chicago IL, USA

Presented at the Annual Meeting of the American Society for Laser Medicine and Surgery, 2018 Dallas, TX.

HIGHLIGHTS

- 16 patients received 5-8 treatments to evaluate effects of an extended protocol. Subject were evaluated 1 month post-treatments.
- Abdominal fat thickness was reduced on average by 19.2 % or 3.4 mm.
- Simultaneously a **15.8** % increase in abdominal muscle thickness was observed, coupled with a 10.8 % reduction in diastatis recti.
- Waist circumference decreased on average by 1.2 inch (after 4th Tx) and 1.6 inch (after the last Tx).
- Data suggest 4 treatments as the ideal protocol.




BASELINE 1 MONTH FU



RESULTS

UMBILICAL CIRCUMFERENCE


• The average circumference decreased by 3.04 cm and 4.17 cm after 4th and last (5th to 8th) treatment, respectively (p<0.003)





CT MEASUREMENTS

CT calculated thickness of rectus abdominis at baseline and 1 month post treatments.

Subcutaneous fat thickness at baseline and 1-month post treatments. Patient ID8 fat measurements could not be objectively made due to close-to-zero baseline fat thickness.







BASELINE 1 MONTH FU





CT scans of patient ID9 at baseline (left) and 1-month post treatments (right). The scan shows reduction of subcutaneous fat (-30.3%) and thickening of rectus abdominis muscle (+8.4%).